INSTRUCTIONS FOR USE

NAME
PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold)

INTENDED USE
The PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) is used for qualitative detection of SARS-CoV-2 virus nucleocapsid protein which is an important conserved structural protein of SARS-CoV-2 in human nasal/oropharyngeal samples. The detection of SARS-CoV-2 nucleocapsid protein antigen can be used to assist the diagnosis of novel coronavirus infection and is helpful for the early detection of novel coronavirus pneumonia infection in the latent period.

TEST PRINCIPLES
The PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (colloidal Gold) is an immune sandwich assay. The test strip contains membranes which are pre-coated with mouse anti-CoV N protein monoclonal antibodies on the test lines. Another mouse anti-CoV N protein monoclonal antibodies which can specifically bind to SARS-CoV-2 N protein are bound to gold particles and sprayed on conjugation pads, when the sample is applied to the sample wells, SARS-CoV N protein and labeled antibody are formed and travel up the strip. The labeled reagent is used to form a visible red line. The presence of SARS-CoV-2 will be indicated by a visible red test line (T) in the result window. Membrane is pre-coated with chicken IgY on the control (C) line. control (c) line appears in each result window when sample has flowed through the strip. The control Line is used as a procedural control. The control line should always appears when the test procedure is performed properly and the reagents are working as intended.

Number	packaging specifications	Test card	Nasal swab Oropharyngeal swab (optional)	Extraction Tube	Test buffer (Extraction Reagent)
1	1 test / box	1	1	1	1
2	5 tests / box	5	5	5	5
3	25 tests / box	25	25	25	25
4	50 tests / box	50	50	50	50

COMPOSITION OF THE REAGENT PRODUCT

1. Reagent device is comprised of a plastic cassette and a test strip. The test strip is comprised of a sample pad, a conjugate pad, nitrocellulose membrane, a wicking pad and polyvinyl chloride(PVC) backing card.
2. Test buffer is composed of PBS buffer, Triton-100 and Proclin 300.
3. Do not exchange product components with different batch numbers.
4. Components not included in the product but necessary for the experiment: Timer.

WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. Do not use the kit contents beyond the expiration date printed on the outside of the box.
3. The test kit has been authorized only for the detection of nucleocapsid proteins from SARS-CoV-2, not for any other viruses or pathogens
4. To avoid erroneous results, specimens must be processed as indicated in the assay procedure section.
5. use appropriate precautions in the collection, handling, storage and disposal of patient samples and used kit contents.
6. Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if operator does not carefully follow the instructions in this package insert.
7. The best detection temperature is 20°C ~ 30°C and humidity RH is 40% ~ 60%.
8. The experiment should be carried out as soon as possible to avoid the inaccurate test result caused by extended exposure of strip in the air after the strip is taken out from the foil pouch.
9. Dispose all specimens and materials resulted from the test as biohazardous waste.
10. Negative results do not rule out SARS-CoV-2 infection, particularly in those who had been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals.
11. Positive results may be due to past or present infection with non- SARS-CoV-2 coronavirus strains.

SPECIMEN COLLECTION AND HGNDLING

Specimen Collection
Anterior Nasal Specimen Collection: Ask the patients to blow their nose. The swab tip should be inserted to a maximum depth of 2.5 cm from the edge of the nostril. Roll the swab over the mucosa of the nostril 5 times to ensure that both mucus and cells are collected. Repeat this procedure with the same swab for the other nostril to ensure that a sufficient amount of sample is collected. Withdraw the swab from the nostril.
Nasopharyngeal Specimen Collection:Ask the patients to blow their nose. Tilt the patient's head back 70 degrees and insert the swab into the nostril parallel to the palate (not upwards) until you encounter resistance or the distance from the nostrils to the outer opening of the patient's ear is equal, indicating contact with the nasopharynx. Leave the swab in place for a few seconds to absorb secretions. Remove the swab slowly, rotating it as you do so.
Oropharyngeal Specim Collection: Insert swab into the posterior pharynx and tonsil areas. This refers to the area around the uvula. Rub the swab over both tonsillar pillars and the posterior oropharynx, avoiding touching the tongue, teeth and gums.

INSTRUCTIONS FOR USE

PERFORMANCE CHARACTERISTICS

1. Clinical Evaluation:
Clinical evaluation was performed to compare the results obtained by novel coronavirus (SARS-CoV-2) Antigen Rapid Test Kit(swab) and PCR. All the observation results and findings in the clinical trial are verified to ensure the reliability of the data and ensure that the conclusions in the clinical trial come from the original data. There are corresponding data management measures in clinical trial and data processing stage. The results were summarized as below:

Consult Instructions For Use		2019-nCoV Nucleic Acid Test Kit (RT-PCR)		Total Results
The Novel coronavirus (SARS-CoV-2) Antigen Rapid Test Kit (Colloidal Gold)	Results	Positive	Negative	
		Positive	98	
	Negative	7	237	244
Total Results		105	239	344

Clinical sensitivity = 98/105=93.33% (95%CI*86.27%~97.05%)
Clinical specificity = 237/239=99.16% (95%CI*96.68%~99.85%)
Accuracy: (98+237)/(98+7+2+237) *100%=97.38% (95%CI* 96.58% to 99.43%)
* Confidence Interval

2. Specimen material: The specimen material for the clinical evaluation is nasopharyngeal swab(NP)and oropharyngeal swab(op)samples.

3. Limit of Detection (LoD): Studies determined the lowest detectable concentration of SARS-CoV-2 at which 100% of all (true positive) replicates test positive. The LoD for the PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal Gold) was established using limiting dilutions of gamma-irradiated SARS-CoV-2. The final LoD was determined to be the lowest concentration resulting in positive detection of twenty(20)out of twenty (20)replicates. Based on this testing the LoD for nasal swab specimens was confirmed as : 121TCID₅₀/mL.

4. High dose hook effect: Strong positive samples of novel coronavirus antigen were diluted in a gradient,then gradient samples from low concentration to high concentration were tested, and were repeated in 5 times. The concentration when the colour depth decreased with increased concentration was used as the minimum concentration of novel coronavirus Antigen when the hook effect appeared. A testing up to 10⁷TCID₅₀/mL of inactivated SARS-CoV-2 to demonstrate that a high dose effect will not occur.

5. Cross-Reactivity: Cross-reactivity of the PocRoc® SARS-CoV-2 Antigen Rapid Test Kit(Colloidal Gold)was evaluated by testing various viruses as below, the final results show that the kit does not cross react with any other viruses and microorganisms except SARS-Cov (2003).

Substances	Concentration	Substances	Concentration
Coronavirus OC43	10 ⁷ PFU/mL	Parainfluenza virus 2	10 ⁷ PFU/mL
Coronavirus NL63	10 ⁷ PFU/mL	Parainfluenza virus 3	10 ⁷ PFU/mL
Coronavirus 229E	10 ⁷ PFU/mL	Parainfluenza virus 4	10 ⁷ PFU/mL
Coronavirus HKU1	10 ⁷ PFU/mL	Human Metapneumovirus (hMPV)	10 ⁷ PFU/mL
Influenza A	10 ⁷ PFU/mL	Adenovirus C1	10 ⁷ PFU/mL
Influenza B	10 ⁷ PFU/mL	Adenovirus 71	10 ⁷ PFU/mL
SARS-coronavirus	10 ⁷ PFU/mL	Rhinovirus A16	10 ⁷ PFU/mL
MERS-coronavirus	10 ⁷ PFU/mL	Streptococcus pneumoniae	10 ⁷ CFU/mL
Respiratory syncytial virus	10 ⁷ PFU/mL	Legionella pneumophila	10 ⁷ CFU/mL
Parainfluenza virus 1	10 ⁷ PFU/mL	Mycoplasma pneumoniae	10 ⁷ CFU/mL


6. Endogenous Interference substances studies: A list of substances that tested for endogenous substances interference study with inactivated SARS-CoV-2 at 3xLoD. The results are shown in the following table ; there was no interference between the detection reagents and the potential interference substances.

Substances	Concentration	Substances	Concentration
Whole Blood	5%(V/V)	Zicam	5%(V/V)


TEST PROCEDURE

1. Open the package and take out the test card.
2. Place the extraction tube pre-filled with test buffer in the holder, then tear the sealing foil of the tube.
3. Remove the swab from its packaging and insert it through the mouth to the wall of the throat or into both nostrils. Follow the instructions for specimen collection.
4. After collecting the specimen, put the swab into the extraction tube, rotate the swab for about 10 seconds and press the swab head against the tube wall to release the antigen in the swab. When removing the swab, squeeze the edge of the tube slightly to remove as much liquid as possible from the tube. Dispose of swabs according to biohazard waste disposal method.
5. Install the plug on the extraction tube, put 3 drops into the sample hole of the test card and start the timer.
6. Read the results after 20 minutes. A strong positive result can be reported within 20 minutes, but a negative result must be reported after 20 minutes and the result after 30 minutes is no longer valid.

INTERPRETATION OF ASSAY RESULT
Interpretation of Results
Negative Results: One color line appears in the control region (c). No apparent red or pink line appears in the test region (T). (As below)

Positive Results: Two distinct colored lines appear. One color line should be in the control region (c) and another color line should be in the test region (T) means positive. (As below)

Invalid Result: If the QC Line C is not observed, the detection should be re-detected regardless of whether or not the detection line is displayed. (As below)

QUALITY CONTROL OF THE PRODUCT

1. Comprehensive internal quality management processes are established to ensure the product quality. The quality of all relevant processes such as of supplier control procedures, production process control procedures, product inspection control procedures, finished product release control procedures is monitored consequently.
2. The measurement results can be traced to the cell culture provided by the supplier. The specific process of quantity value transfer is reflected in the traceability data in the technical documents.

LIMITATIONS

1. The PocRoc® SARS-CoV-2 Antigen Rapid Test Kit (Colloidal gold) does not differentiate between SARS-CoV and SARS-CoV-2.
2. This reagent is only used for qualitative testing and does not indicate the number of novel coronavirus antigen in the sample.
3. This reagent is only a clinical auxiliary diagnostic tool. If the results are positive, it is recommended to use other methods for further examination and to follow the doctor's diagnosis.
4. The test results of this product are for clinical reference only and should not be used as the sole basis for clinical diagnosis and treatment. The clinical management of patients should be combined with their symptoms, signs, medical history, other laboratory examinations, treatment response and epidemiology.
5. In case of changes of the analytical performance of the product, the customer will be timely informed according to the customer notification procedure of the manufacturer.

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